

IT'S MORE THAN
3D PRINTING.



NEW

IT'S
HIGHER PRODUCTIVITY.

CARA PRINT CUBE

Want to simplify your workflow and improve patient care?

The cara Print Cube is your compact, reliable solution for seamless 3D printing. It delivers high precision with an expanded build area for up to five dental models, ensuring consistent, exceptional results.

Cost-effective and easy to use, this 3D printer is designed for dental professionals and fits effortlessly into labs, practices and educational institutions.

Experience smarter dental manufacturing with the cara Print Cube – where simplicity and productivity meet.



kulzer.com/mycube



KULZER
MITSUI CHEMICALS GROUP

CARA PRINT CUBE

- ▶ **Expanded build area:** Tailored for 5 dental models
- ▶ **High precision:** 6K resolution and 34 µm accuracy.
- ▶ **Cost-Effective:** Competitive price with superior performance.
- ▶ **Validated system:** Compatible with dima Print materials, ensuring seamless integration.



Main parameters

Build platform/volume (l x w x h) [mm]	192 x 120 x 190
Pixel size [µm]	34.4
Technology	6K LCD
Layer thickness [mm]	0.05 – 0.3
Light Source	405 nm LED
Resolution	5760 x 3600 pixel
Input voltage	100 ~ 120 VAC, 4 A, 50/60 Hz 200 ~ 240 VAC, 2 A, 50/60 Hz
Device Dimensions (l x w x h) [mm]	390 x 420 x 550
Device Net Weight [kg]	31

