

DIMA PRINT PERFORMANCE OVERVIEW


| PICTURE | MATERIAL | INDICATION | EU MD CLASS | COLOUR | PRINTING TIME | EQUIVALENT PRINTING TIME PER PART * | COST PER PART ** |
|---|--------------------------|--|----------------|--|---|-------------------------------------|------------------|
|  | dima Print C&B temp | Temporary crowns & bridges | Ila | A1, A2, A3, A3.5, B1, B2, BL1 | 1 h 34 min (for 110 units) | 3 min | 0,38 € |
|  | dima Print Cast ruby | CAD to cast | not applicable | red | 1 h 28 min (for various CAD-to-Cast objects) | 5 min | 0,28 € |
|  | dima Print Denture Base | Full and partial removable dentures and baseplates | Ila | original pink, light pink, light reddish pink, dark pink | 1 h 44 min (for 10 arches) | 12 min | 5,91 € |
|  | dima Print Denture Teeth | Denture Teeth | Ila | A1, A2, A3, A3.5, B1, B2, BL1 | 1 h 35 min (for 8 arches) | 4 min | 0,28 € |
|  | dima Print Splint clear | Splint / Nightguards | Ila | transparent clear | 1 h 18 min (for 6 Splints) | 13 min | 2,25 € |
|  | dima Print Stone beige | Dental Models | not applicable | gypsum beige | 1 h 30 min (for 4 big models) | 22 min | 12,47 € |

Availability in your country may be restricted due to legal requirements. Please check the availability for your country with Kulzer Customer Service Center.

* depending on various factors, especially on the utilisation of the build platform; calculation based on build platform utilisation as shown in picture

** based on the RRP of the respective dima Print material only

Find the current version here: kulzer.com/mycube-support



KULZER
MITSUI CHEMICALS GROUP

DIMA PRINT PERFORMANCE OVERVIEW

| PICTURE | MATERIAL | INDICATION | EU MD CLASS | COLOUR | PRINTING TIME | EQUIVALENT PRINTING TIME PER PART * | COST PER PART ** |
|---|--------------------------------|------------------------------------|----------------|--------------------------|-------------------------------------|-------------------------------------|------------------|
|  | dima Print Stone beige | Dental Models | not applicable | gypsum beige | 1 h 18 min (for 5 models) | 13 min | 9,54 € |
|  | dima Print Stone beige | Dental Models | not applicable | gypsum beige | 1 h 15 min (for 1 model) | 1 h 15 min | 9,54 € |
|  | dima Print Guide & Tray | Impression Trays & Surgical Guides | I | transparent light orange | 1 h 43 min (for 5 trays) | 28 min | 6,25 € |
|  | dima Print Guide & Tray | Impression Trays & Surgical Guides | I | transparent light orange | 57 min (for 6 big guides) | 10 min | 2,32 € |
|  | dima Print Denture Base Try-In | Try-ins | IIa | white | 1 h 01 min (for 6 arches) | 10 min | 9,79 € |
|  | dima Print Stone beige Smart | Dental Models | not applicable | gypsum beige | 21 min (for 1 model for splints) | 21 min | 3,53 € |

Availability in your country may be restricted due to legal requirements. Please check the availability for your country with Kulzer Customer Service Center.

* depending on various factors, especially on the utilisation of the build platform; calculation based on build platform utilisation as shown in picture

** based on the RRP of the respective dima Print material only

Find the current version here: kulzer.com/mycube-support



DIMA PRINT PERFORMANCE OVERVIEW

| PICTURE | MATERIAL | INDICATION | EU MD CLASS | COLOUR | PRINTING TIME | EQUIVALENT PRINTING TIME PER PART * | COST PER PART ** |
|---|------------------------------|--|----------------|-------------------------------|--------------------------------------|-------------------------------------|------------------|
|  | dima Print Stone beige Smart | Dental Models | not applicable | gypsum beige | 23 min (for 6 models for splints) | 4 min | 3,53 € |
|  | dima Print Denture Base 2 | Full and partial removable dentures and baseplates | IIa | orange-pink | 1h 54 min (for 10 arches) | 13 min | 5,07 € |
|  | dima Print Gingiva Mask | Gingiva Masks | not applicable | pink | 17 min (for 10 gingiva masks) | 2 min | 1,42 € |
|  | dima Print Stone teal | Dental models for deep drawing | not applicable | turquoise | 1h 57 min (for 5 models) | 23 min | 10,30 € |
|  | dima Print Denture Teeth 2 | Temporary Crowns & Bridges | IIa | A1, A2, A3, A3.5, B1, B2, BL1 | 1h 59 min | 4 min | 0,93 € |
|  | dima Print Denture Teeth 2 | Denture Teeth | IIa | A1, A2, A3, A3.5, B1, B2, BL1 | 1h 56 min | 5 min | 1,12 € |

Availability in your country may be restricted due to legal requirements. Please check the availability for your country with Kulzer Customer Service Center.

* depending on various factors, especially on the utilisation of the build platform; calculation based on build platform utilisation as shown in picture

** based on the RRP of the respective dima Print material only

Find the current version here: kulzer.com/mycube-support

