







## DIMA PRINT PERFORMANCE OVERVIEW

PICTURE	MATERIAL	INDICATION	MD CLASS	COLOUR	PRINTING TIME	EQUIVALENT PRINTING TIME PER PART*	COSTS PER PART [€]**
	dima Print C&B temp	Temporary Crowns & Bridges, total of 102 units	Ila	A1, A2, A3, A3.5, B1, B2, BL1	1 h 34 min	3 min	0,38€
	dima Print Cast ruby	CAD to cast, total of 19 objects	not applicable	red	1 h 28 min	5 min	0,28€
	dima Print Denture Base	Full denture bases, 10 arches	Ila	original pink, light pink, light reddisch pink, dark pink	1 h 44 min	12 min	5,91€
	dima Print Denture Teeth	Denture Teeth	Ila	A1, A2, A3, A3.5, B1, B2, BL1	1 h 35 min	4 min	0,28€
	dima Print Splint clear	6 Splints/ Nightguards	Ila	transparent clear	1 h 18 min	13 min	2,25€
	dima Print Stone beige	4 big Dental models with removable dies	not applicable	gypsum beige	50 µm: 3 h 27 min 100 µm: 1 h 30 min	50 µm: 52 min 100 µm: 22 min	12,47€

Availability in your country may be restricted due to legal requirements. Please check the availability for your country with Kulzer Customer Service Center.

\*depending on various factors, especially on the utilization of the build platform, calculation based on build platform utilization as shown in picture

\*\*based on RRP

Current version find here: [kulzer.com/mycube-support](https://kulzer.com/mycube-support)



**KULZER**  
MITSUI CHEMICALS GROUP

## DIMA PRINT PERFORMANCE OVERVIEW

PICTURE	MATERIAL	INDICATION	MD CLASS	COLOUR	PRINTING TIME	EQUIVALENT PRINTING TIME PER PART*	COSTS PER PART [€]**
	dima Print Stone beige	5 models for splint	not applicable	gypsum beige	50 µm: 2 h 52 min 100 µm: 1 h 15 min	50 µm: 34 min 100 µm: 15 min	9,54€
	dima Print Stone beige	one model for splint	not applicable	gypsum beige	50 µm: 3 h 9 min 100 µm: 1 h 18 min	50 µm: 3 h 9 min 100 µm: 1 h 18 min	10,32€
	dima Print Guide & Tray	5 Impression trays	I	transparent light orange	150 µm: 1 h 43 min 100 µm: 2 h 20 min	150 µm: 21 min 100 µm: 28 min	7,15€
	dima Print Try-in	Try-ins, 6 arches	Ila	white	1 h 01 min	10 min	9,79€
	dima Print Stone beige Smart	one model for splint	not applicable	gypsum beige	100 µm: 14 min	14 min	3,53€
	dima Print Stone beige Smart	6 models for splint	not applicable	gypsum beige	100 µm: 16 min	3 min	3,53€

Availability in your country may be restricted due to legal requirements. Please check the availability for your country with Kulzer Customer Service Center.

\*depending on various factors, especially on the utilization of the build platform, calculation based on build platform utilization as shown in picture

\*\*based on RRP

Current version find here: [kulzer.com/mycube-support](https://kulzer.com/mycube-support)

