

## CARA PRINT CUBE - WEAR PARTS



ARTICLE CODE	NAME	DISCRIPTION	CATEGORY	PRODUCT IMAGE
66099965	Resin tank; Cube	Resin tank incl. cover	Wear part	
66099975	Build platform; Cube	Build platform for cara Print Cube	Wear part	
66099976	LCD screen; Cube	LCD screen incl. connection cable	Wear part	
66099971	Release Film Kit; Cube	Spare release films for the resin tank incl. accessories, screws, tapes and LCD screen protection foils	Wear part	





# CARA PRINT CUBE - WEAR PARTS

ARTICLE CODE	NAME	DISCRIPTION	CATEGORY	PRODUCT IMAGE
66099972	Air filter; Cube	Spare air filter for the cara Print Cube printer	Wear part	
66099973	Tool Kit; Cube	Incl. All relevant helping tools for operation and maintenance	Wear part	
66099974	Tray kit; Cube	Incl. tray and hammer	Wear part	
66099977	Cover hood; Cube	Spare cover hood	Wear part	





# CARA PRINT CUBE - WEAR PARTS

ARTICLE CODE	NAME	DISCRIPTION	CATEGORY	PRODUCT IMAGE
66100104	Resin tank fastening; Cube	Sspare resin tank fastening incl. screws. Always comes a a pair of two pieces.	Wear part	
66099978	Fresnel lens; Cube	Spare fresnel lens, below LCD screen, below glass substrate	Wear part	
66100107	Glass substrate; Cube	Spare glass substrate, below LCD screen, above fresnel lens	Wear part	
66100149	CoExtra license; Cockpit	Additional life-long cara Print Cockpit license, virtual article	Software	



Find the current version here: [kulzer.com/mycube-support](https://kulzer.com/mycube-support)



## CARA PRINT CLEANUP - WEAR PARTS

ARTICLE CODE	NAME	DISCRIPTION	CATEGORY	PRODUCT IMAGE
66100120	Washing bucket; CleanUp	Washing bucket for cara Print CleanUp, w/o metal basket	Wear part	





## CARA PRINT LEDCURE - WEAR PARTS

ARTICLE CODE	NAME	DISCRIPTION	CATEGORY	PRODUCT IMAGE
66094927	Turntable, transparent; LED	Transparent turntable	Wear part	

