

## artic digital

GB Mould chart  
DE Formenkarte  
FR Carte des Formes  
ES Carta de Formas  
IT Tabella forme  
PT Carta Molde  
NL Vormenkaart  
SE Formkort  
DK Formkort  
FI Muotto korti  
NO Formkort  
EE Vormikaart  
LV Formas karte  
LT Formos kortelė  
PL Karta kształtu  
IS Myglatöflu  
HR Kartica oblika  
BG плесен диаграма  
CZ Karta tvaru  
SK Tvarová karta  
RS картица калупа  
MK форма картичка  
KR 모양 카드  
GR Πίνακας σχημάτων  
ARA مخطط القالب



**Artic® Digital**

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP

## artic digital 6

Anteriors / Frontzähne / Antérieures / Anteriores / Anteriori / Anteriores / Tanden /  
 Framtänder / For tænder / Etuhampaat / Fortennene / Eesmised hambad /  
 Priekšējie zobi / Priekiniai dantys / Zęby przednie / Framhliðar / Prednji zubi / предни зъби /  
 Přední zuby / Predné zuby / предњи зуби / предните заби / 앞니 / Πρόσθια / الأمامية

|                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                           |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>O 20</b><br>   <div> <div>46,1</div> <div>9,7</div> <div>8,4</div> </div>        | <b>S 35</b><br>   <div> <div>50,6</div> <div>10,8</div> <div>9,2</div> </div>       |
| <b>O 30</b><br>   <div> <div>46,6</div> <div>10,9</div> <div>8,4</div> </div>       | <b>S 40</b><br>   <div> <div>52,3</div> <div>12,4</div> <div>9,8</div> </div>       |
| <b>R 20</b><br>   <div> <div>42,9</div> <div>9,9</div> <div>7,8</div> </div>        | <b>T 10</b><br>   <div> <div>44,1</div> <div>9,0</div> <div>8,2</div> </div>        |
| <b>R 30</b><br>   <div> <div>47,7</div> <div>12,1</div> <div>8,8</div> </div>   | <b>T 20</b><br>   <div> <div>45,6</div> <div>9,8</div> <div>8,8</div> </div>    |
| <b>R 40</b><br>   <div> <div>51,1</div> <div>13,3</div> <div>9,5</div> </div> | <b>T 30</b><br>   <div> <div>45,0</div> <div>10,8</div> <div>8,1</div> </div> |
| <b>S 05</b><br>   <div> <div>40,3</div> <div>8,9</div> <div>7,4</div> </div>  | <b>T 35</b><br>   <div> <div>47,4</div> <div>10,0</div> <div>8,8</div> </div> |
| <b>S 10</b><br>   <div> <div>44,2</div> <div>9,1</div> <div>7,9</div> </div>  | <b>T 40</b><br>   <div> <div>49,1</div> <div>12,1</div> <div>8,7</div> </div> |
| <b>S 20</b><br>   <div> <div>44,0</div> <div>9,9</div> <div>7,9</div> </div>  | <b>TO 05</b><br>   <div> <div>42,1</div> <div>9,0</div> <div>7,8</div> </div> |
| <b>S 30</b><br>   <div> <div>48,2</div> <div>10,9</div> <div>8,6</div> </div> | <b>TO 10</b><br>   <div> <div>43,3</div> <div>9,3</div> <div>8,3</div> </div> |


## artic digital 6

Anteriors / Frontzähne / Antérieures / Anteriores / Anteriori / Anteriores / Tanden /  
 Framtänder / For tænder / Etuhampaat / Fortennene / Eesmised hambad / Priekšējie  
 zobi / Priekiniai dantys / Zęby przednie / Framhliðar / Prednji zubi / предни зъби /  
 Přední zuby / Predné zuby / предњи зуби / предните заби / 앞니 / Πρόσθια / الأمامية


|                                                                                                                                                                                                                                                                       |                                                |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------|
| <b>TO 20</b><br>   | <div>43,6</div> <div>10,2</div> <div>7,7</div> |
| <b>L 25</b><br>                                                                                  | <div>36,8</div> <div>9,3</div> <div>5,5</div>  |
| <b>TO 30</b><br>   | <div>45,9</div> <div>10,4</div> <div>8,6</div> |
| <b>L 30</b><br>                                                                                  | <div>35,9</div> <div>9,0</div> <div>5,5</div>  |
| <b>TO 40</b><br>   | <div>47,8</div> <div>11,0</div> <div>8,9</div> |
| <b>L 40</b><br>                                                                                  | <div>36,0</div> <div>9,0</div> <div>5,5</div>  |
| <b>L 10</b><br>                                                                                   | <div>33,5</div> <div>7,8</div> <div>5,1</div>  |
| <b>L 45</b><br>                                                                                | <div>38,9</div> <div>9,6</div> <div>5,9</div>  |
| <b>L 15</b><br>                                                                                 | <div>33,9</div> <div>9,5</div> <div>5,2</div>  |
| <b>L 50</b><br>                                                                              | <div>39,2</div> <div>11,1</div> <div>5,9</div> |
| <b>L 20</b><br>                                                                                 | <div>35,2</div> <div>9,1</div> <div>5,4</div>  |
| <b>L 55</b><br>                                                                              | <div>40,2</div> <div>10,6</div> <div>6,0</div> |

## artic digital 8 0°/20°

Posteriors/Seitenzähne/Postérieures/Posteriores/Posteriori/Posteriores/Kiezen/Bakre  
tänder/Bageste tænder/Takimmaisest hampaat/Bakre tenner/Tagumised hambad/  
Aizmugurējīe zobi/Užpakalīnīai dantys/Zęby tylne/Bakhliðar/Stražņji zubi/задни  
зъби/Zadní zuby/Zadné zuby/задњи зуби/задните заби/구치부/Οπίσθια/الخلفية

|      |     |                                                                                      |                                                                                       |      |
|------|-----|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|------|
| 30 U | 0°  |     |    | 27,2 |
| 30 L |     |                                                                                      |                                                                                       | 7,9  |
| 32 U | 0°  |    |    | 33,1 |
| 32 L |     |                                                                                      |                                                                                       | 8,8  |
| 32 U | 20° |  |  | 31,1 |
| 32 L |     |                                                                                      |                                                                                       | 9,1  |
| 34 U | 20° |  |  | 32,9 |
| 34 L |     |                                                                                      |                                                                                       | 9,4  |
| 36 U | 20° |  |  | 34,7 |
| 36 L |     |                                                                                      |                                                                                       | 10,0 |
|      |     |                                                                                      |                                                                                       | 37,4 |
|      |     |                                                                                      |                                                                                       | 9,0  |

|                                                                                     | Upper anteriors                                                                           | Lower anteriors | Posteriors | 0°     | 20°    |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|-----------------|------------|--------|--------|
|    |  O 20    | L 30            |            | 30, 32 | 32     |
|                                                                                     |  O 30    | L 40            |            | 32     | 32     |
|    |  R 20    | L 20            |            | 32     | 32, 34 |
|                                                                                     |  R 30    | L 25            |            | 32     | 36     |
|                                                                                     |  R 40    | L 55            |            | 32     | 36     |
|    |  S 05    | L 10            |            | 30     | 32     |
|                                                                                     |  S 10    | L 10            |            | 30     | 32     |
|                                                                                     |  S 20    | L 20            |            | 30, 32 | 32     |
|                                                                                     |  S 30    | L 40            |            | 32     | 32, 34 |
|                                                                                     |  S 35    | L 55            |            | 32     | 36     |
|                                                                                     |  S 40    | L 55            |            | 32     | 36     |
|                                                                                     |                                                                                           |                 |            |        |        |
|    |  T 10    | L 10            |            | 30     | 32     |
|                                                                                     |  T 20    | L 30            |            | 30, 32 | 32     |
|                                                                                     |  T 30    | L 20            |            | 32     | 32     |
|                                                                                     |  T 35    | L 40            |            | 32     | 32     |
|                                                                                     |  T 40   | L 50            |            | 32     | 34     |
|  |  TO 05 | L 10            |            | 30     | 32     |
|                                                                                     |  TO 10 | L 15            |            | 30     | 32     |
|                                                                                     |  TO 20 | L 20            |            | 30, 32 | 32     |
|                                                                                     |  TO 30 | L 40            |            | 32     | 32     |
|                                                                                     |  TO 40 | L 45            |            | 32     | 34     |



This mould chart complies with EN ISO 22112:2017. / Diese Formenkarte erfüllt EN ISO 22112:2017. / Cette carte forme remplit EN ISO 22112:2017. / Esta carta de formas es conforme con la EN ISO 22112:2017. / Questa tabella forme soddisfa EN ISO 22112:2017. / Esta carta molde está em conformidade com EN ISO 22112:2017. / Deze vormenkaart voldoet aan EN ISO 22112:2017. / Detta formkort uppfyller EN ISO 22112:2017. / Dette formkort opfylder EN ISO 22112:2017. / Tämä muotokortti täyttää EN ISO 22112:2017. / Dette formkortet oppfyller EN ISO 22112:2017. / See vormi-kaart vastab standardile EN ISO 22112:2017. / Ši formas karte atbilst EN ISO 22112:2017. / Ši formos korta išpildo EN ISO 22112:2017. / Ta karta kształtu spełnia EN ISO 22112:2017. / Petta mótakort er í samræmi við EN ISO 22112:2017. / Ova kartica oblika ispunjava EN ISO 22112:2017. / Тази матрица е в съответствие с EN ISO 22112:2017. / Tato karta tvaru splňuje EN ISO 22112:2017. / Táto tvarová karta spĺňa EN ISO 22112:2017. / Ова картица калупа је у складу са ЕН ИСО 22112:2017. / Оваа карта за мувла е во согласност со EN ISO 22112:2017. / 이 몰드 카드는 EN ISO 22112:2017을 준수합니다. / Αυτός ο πίνακας σχημάτων συμφωνεί με το πρότυπο EN ISO 22112:2017. / يتوافق مخطط القوالب هذا مع EN ISO 22112:2017.



#### Manufacturer

Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau  
Germany

#### Contact in Germany

Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau  
Germany  
info.lab@kulzer-dental.com

