

artic digital

GB Mould chart
DE Formenkarte
FR Carte des Formes
ES Carta de Formas
IT Tabella forme
PT Carta Molde
NL Vormenkaart
SE Formkort
DK Formkort
FI Muotto korti
NO Formkort
EE Vormikaart
LV Formas karte
LT Formos kortelė
PL Karta kształtu
IS Myglatöflu
HR Kartica oblika
BG плесен диаграма
CZ Karta tvaru
SK Tvarová karta
RS картица калупа
MK форма картичка
KR 모양 카드
GR Πίνακας σχημάτων
ARA مخطط القالب



Artic® Digital









































Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

artic digital 6

Anteriors/ Frontzähne/ Antérieures/ Anteriores/ Anteriori/ Anteriores/ Tanden/ Framtänder/ For tænder/ Etuhampaat/ Fortennene/ Eesmised hambad/ Priekšējie zobi/ Priekiniai dantys/ Zęby przednie/ Framhliðar/ Prednji zubi/ предни зъби/ Přední zuby/ Predné zuby/ предњи зуби/ предните зуби/ 앞니/ Πρόσθια/ الأمامية

O10D  	40,9 9,0 7,5	S35D  	50,6 10,8 9,2
O20D  	46,1 9,7 8,4	S40D  	52,3 12,4 9,8
O30D  	46,6 10,9 8,4	T10D  	44,1 9,0 8,2
R20D  	42,9 9,9 7,8	T15D  	42,6 9,6 7,9
R30D  	47,7 12,1 8,8	T20D  	45,6 9,8 8,8
R40D  	51,1 13,3 9,5	T30D  	45,0 10,8 8,1
S05D  	40,3 8,9 7,4	T35D  	47,4 10,0 8,8
S10D  	44,2 9,1 7,9	T40D  	49,1 12,1 8,7
S20D  	44,0 9,9 7,9	TO05D  	42,1 9,0 7,8
S30D  	48,2 10,9 8,6	TO10D  	43,3 9,3 8,3

artic digital 6

Anteriors/ Frontzähne/ Antérieures/ Anteriores/ Anteriori/ Anteriores/ Tanden/ Framtänder/ For tænder/ Etuhampaat/ Fortennene/ Eesmised hambad/ Priekšējie zobi/ Priekiniai dantys/ Zęby przednie/ Framhliðar/ Prednji zubi/ предни зъби/ Přední zuby/ Predné zuby/ предњи зъби/ предните заби/ 앞니/ Πρόσθια/ الأمامية














artic digital 8 0°

Posteriors/ Seitenzähne/ Postérieures/ Posteriores/ Posteriori/ Posteriores/ Kiezen/ Bakre tänder/ Bageste tænder/ Takimmaisest hampaat/ Bakre tenner/ Tagumised hambad/ Aizmugurējie zobi/ Užpakaliniai dantys/ Zęby tylne/ Bakhliðar/ Stražnji zubi/ задни зъби/ Zadní zuby/ Zadné zuby/ задњи зъби/ задните заби/ 구치부/ Οπίσθια/ الخلفية


TO 20 D		43,6 10,2 7,7	L 25 D		36,8 9,3 5,5
TO 30 D		45,9 10,4 8,6	L 30 D		35,9 9,0 5,5
TO 40 D		47,8 11,0 8,9	L 40 D		36,0 9,0 5,5
L 05 D		31,5 9,0 5,0	L 45 D		38,9 9,6 5,9
L 10 D		33,5 7,8 5,1	L 50 D		39,2 11,1 5,9
L 15 D		33,9 9,5 5,2	L 55 D		40,2 10,6 6,0
L 20 D		35,2 9,1 5,4			
30 UD	0°				27,2 7,9
30 LD					28,9 8,6
32 UD	0°				33,1 8,8
32 LD					34,8 9,2

artic digital 8 10°/20°

Posteriors/Seitenzähne/Postérieures/Posteriores/Posteriori/Posteriores/Kiezen/Bakre tänder/Bageste tænder/
 Takimmaisest hampaat/Bakre tenner/Tagumised hambad/Aizmugurējie zobi/Užpakaliniai dantys/Zęby tylne/Bakhliðar/
 Stražnji zubi/задни зъби/Zadní zuby/Zadné zuby/задњи зуби/задните заби/구치부/Οπίσθια/الخلفية

30 UD	10°			29,4
30 LD				8,8
32 UD	10°			32,1
32 LD				9,5
30 UD	20°			29,2
30 LD				7,6
32 UD	20°			31,0
32 LD				6,9
34 UD	20°			31,1
34 LD				9,1
36 UD	20°			33,9
36 LD				8,0
34 UD	20°			32,9
34 LD				9,4
36 UD	20°			35,7
36 LD				8,3
36 UD	20°			34,7
36 LD				10,0
36 UD	20°			37,4
36 LD				9,0

	Upper anteriors	Lower anteriors	Posteriors	0°	10°	20°
○	● O 10 D	L 10 D		30	30	30
	● O 20 D	L 30 D		30, 32	30, 32	32
	● O 30 D	L 40 D		32	32	32
◻	■ R 20 D	L 20 D		32	32	32, 34
	■ R 30 D	L 25 D		32	32	36
	■ R 40 D	L 55 D		32	32	36
◻	■ S 05 D	L 10 D		30	30	30, 32
	■ S 10 D	L 10 D		30	30	30, 32
	■ S 20 D	L 20 D		30, 32	30, 32	32
	■ S 30 D	L 40 D		32	32	32, 34
	■ S 35 D	L 55 D		32	32	36
	■ S 40 D	L 55 D		32	32	36
△	▲ T 10 D	L 10 D		30	30	30, 32
	▲ T 15 D	L 15 D		30, 32	30, 32	32
	▲ T 20 D	L 30 D		30, 32	30, 32	32
	▲ T 30 D	L 20 D		32	32	32
	▲ T 35 D	L 40 D		32	32	32
	▲ T 40 D	L 50 D		32	32	34
◻	▲ TO 05 D	L 05 D		30	30	30
	▲ TO 05 D	L 10 D		30	30	30, 32
	▲ TO 10 D	L 15 D		30	30	30, 32
	▲ TO 20 D	L 20 D		30, 32	30, 32	32
	▲ TO 30 D	L 40 D		32	32	32
	▲ TO 40 D	L 45 D		32	32	34



This mould chart complies with EN ISO 22112:2017. / Diese Formenkarte erfüllt EN ISO 22112:2017. / Cette carte forme remplit EN ISO 22112:2017. / Esta carta de formas es conforme con la EN ISO 22112:2017. / Questa tabella forme soddisfa EN ISO 22112:2017. / Esta carta molde está em conformidade com EN ISO 22112:2017. / Deze vormenkaart voldoet aan EN ISO 22112:2017. / Detta formkort uppfyller EN ISO 22112:2017. / Dette formkort opfylder EN ISO 22112:2017. / Tämä muotokortti täyttää EN ISO 22112:2017. / Dette formkortet oppfyller EN ISO 22112:2017. / See vormi-kaart vastab standardile EN ISO 22112:2017. / Ši formas karte atbilst EN ISO 22112:2017. / Ši formos korta išpildo EN ISO 22112:2017. / Ta karta kształtu spełnia EN ISO 22112:2017. / Þetta mótakort er í samræmi við EN ISO 22112:2017. / Ova kartica oblika ispunjava EN ISO 22112:2017. / Тази матрица е в съответствие с EN ISO 22112:2017. / Tato karta tvaru splňuje EN ISO 22112:2017. / Táto tvarová karta spĺňa EN ISO 22112:2017. / Ова картица калупа је у складу са ЕН ИСО 22112:2017. / Оваа карта за мувла е во согласност со EN ISO 22112:2017. / 이 몰드 카드는 EN ISO 22112:2017을 준수합니다. / Αυτός ο πίνακας σχημάτων συμφωνεί με το πρότυπο EN ISO 22112:2017. / يتوافق مخطط القوالب هذا مع EN ISO 22112:2017.



Manufacturer

Kulzer GmbH
Leipziger Straße 2
63450 Hanau
Germany

Contact in Germany

Kulzer GmbH
Leipziger Straße 2
63450 Hanau
Germany
info.lab@kulzer-dental.com

